HIGG medium

References

Original recipe:

Cohen-Bazire, G., Sistrom, W. R., Stanier, R. Y. 1957. Kinetic studies of pigment by non-sulfur purple bacteria. J. Cell Comp. Physiol. 49:25-68.

Modified for Caulobacter:

Poindexter, J. S. 1978. Selection for nonbuoyant morphological mutants of Caulobacter crescentus. J. Bact. 135:1141-1145.

Phosphate Buffer:

Stoll, V. S., Blanchard, J. S. 1990. Buffers: Principles and Practice. from Guide to Protein Purification. Murray P. Deutscher, ed. Academic Press, Inc. San Diego.

Stock solutions needed

Metals "44" per 100 ml (for addition to Hutner base) 2 1.3 ml O.TM EDTA Disodium EDTA 250.0 mg ZnSO4·7H2O 1095.0 mg FeSO4·7H2O 500.0 mg MnSO₄·H₂O 154.0 mg CuSO4.5H2O 39.2 mg 24.8 Co(NO3)2.6H2O 17.7 mg Na2B4O7.10H2O A few drops of sulphuric acid are added to retard precipitation.

Concentrated Hutner Base (100%) (referred to in Poindextor paper as "double strength") (* omit vitamins if using for Caulobacter):

Nitrilotriacetic acid MgSO₄ CaCl₂·2H₂O (NH4)6M07O24·4H2O FeSO4·7H2O *Nicotinic acid *Thiamin HCl *Biotin Metals "44" Distilled H₂O to 1000 ml Sterile filter (do not autoclave)

5.6 9

KOH

10.0 g herta hydrak anhydrous 595 14.45 g --> 29 72.518 S 3.335 g -9.25 mg 99.0 mg

When making concentrated base, first disolve the nitrilotriacetic acid and neutralize with KOH (it takes approximately 10 ml of 10N solution). Add the rest of the ingredients and adjust pH to 6.6-6.8 before adjusting to volume. This is very important or your Fe will precipitate out. (You will see some precipitation as you adjust the pH but things should go back in solution).

Na2HPO4-KH2PO4 buffer, pH 7.0

0.5 M Phosphate Buffer, pH= 7.0 Be careful to use the correct amount of chemical as these come in multiple forms (anhydrous or hydrated) Calculate the amount of dry chemical you will need to make 0.5 M solutions from the molecular weight of the form of phosphate which you have to use)

50 mg

25.0 mg 0.5 mg

50.0 ml

-61 ml 0.5 M Na₂HPO₄ -39 ml 0.5 M KH2PO4 The pH should be 7.0 upon adding these two components together to make 100 ml; Autoclave.

1 M imidazole pH=7.0 (adjusted with HCl)

20 % glucose (filter sterilized)

20 % sodium glutamate at pH=7.0 (filter sterilized)

1 M NH₄Cl, autoclaved

HIGG medium (1 liter):

The most convenient way to make HIGG is to autoclave the required amount of water, let it cool, and then add each component.

Stock	Volume added for 100 ml	Volume added for 250 ml	Volume added for 500 ml	Volume added for 1 liter	Final concentration
1M imidazole, pH=7	0.5 ml	1.25 ml	2.5 ml	5 ml	5 mM
Hutner concentrated base (100 %)	2 ml	5 ml	10 ml	20 ml	2 %
50mM CaCl2	2 ml	5 ml	10 ml	20ml	1mM
20% glucose	0.75 ml	1.88 ml	3.75 ml	7.5 ml	0.15 %
Nitrogen sources (use both)					Super Constant and
20 % sodium glutamate	0.75 ml	1.88 ml	3.75 ml	7.5 ml	0.15 % (8.9 mM)
1 M NH ₄ Cl	0.89 ml	2.23 ml	4.45 ml	8.9 ml	8.9 mM
Phosphate (500 mM)	6 µL	15 μL	30 µL	60 µL	30 µM (low)
(choose one)	40 µL	100 μL	200 µL	400 µL	0.2 mM
		50 50			balanced P and C
	200 µL	0.5 ml	1 ml	2 ml	1 mM (high)
	2 ml	5 ml	10 ml	20 ml	10 mM (high)

* the published HIGG recipe uses sodium glutamate as a nitrogen source; however, according to Jeanne Poindextor the recipe should include equimolar amounts of both sodium glutamate and NH₄Cl

Doubling time when phosphate is in excess should be approx. 2.5 hours.

Phosphate becomes growth yield limiting at 100 μ M and growth rate limiting at 20 μ M in the presence of 0.2 % glucose. 1 mM phosphate is excess phosphate (non-limiting) and 10 mM can also be used.

When cells are limited for phosphate with respect to carbon, they accumulate poly-β-hydroxybutiric acid (PHB) granules; when cells are limited for carbon with respect to phosphate, they accumulate polyphosphate. Simultaneous exhaustion of phosphate and carbon result in cells lacking either storage polymer [this occurs at 0.2 mM phosphate and 0.2 mM (glucose + glutamate)].